



PrimeSafe SteriliT

Chemical sterilant for the sterilisation of pre-cleaned instruments



Available sizes:
1L; 5L; 20L
Dilution Factor:
Ready to Use

Reg. No: Act5GNR529/242613/110/0601

Locally manufactured: No supply delays

Economical: 7 day usage window, confirmable with MEC Test Strips.

Easy: Ready to use, no extra fuss and no 'cooling down' time.

Protective: Corrosion Inhibitors protect instruments even when not immersed.

Please note:

Where access to an approved, regularly maintained and tested autoclave is available, the autoclave should always be considered the primary means of sterilisation. **PrimeSafe SteriliT** should then be used as a holding medium to maintain the chain of sterility, as well as to ensure the sterility of non-autoclavable and or heat-sensitive items.

For use by healthcare professionals in medical device reprocessing. Recommended for surgical instruments (both medical and dental), endoscopes and any other medical devices which require sterilisation, particularly those that are heat-sensitive.

Antimicrobial Efficacy:

PrimeSafe SteriliT has been tested by the SABS, according to SANS1593:1994 and found to comply with all specifications of the standards.



Storage Instructions:

Store between 5°C & 50°C.
Store out of reach of direct sunlight.
Keep locked out of reach of children at all times.
Store away from foodstuffs.
Consult Material Safety Data Sheet.

Active Ingredient:

PrimeSafe SteriliT contains Gluteraldehyde at a level of 2.75% & QAC at a level of 0.6%.

First Aid:

In case of eye contact:

Flush eye with cool running water for 10-15 minutes. Seek medical advice if irritation persists longer than 30 minutes.

In case of skin contact:

Wash affected area thoroughly with soap and water.

In case of ingestion: Do not induce vomiting, drink a large glass of cold milk and seek medical advice if nausea persists longer than 30 minutes.

In case of inhalation:

Remove affected person to fresh air for at least 10 minutes and seek medical advice if shortness of breath occurs.

Instructions for Use:

Decant the required amount of undiluted PrimeSafe SteriliT into a clean, dry, covered receptacle. Completely immerse thoroughly pre-cleaned, rinsed & dried instruments in undiluted PrimeSafe SteriliT for no less than four (4) hours. Remove the instruments from the liquid using the established post-sterilisation procedure of the institution, and rinse in sterile water prior to re-use. Replace the liquid in the receptacle every 7 days, unless the daily use of M.E.C. Strips indicates sufficient remaining active ingredient levels.

Recommended Uses:

For the chemical sterilisation of pre-cleaned medical devices.

PPE Requirements:

Ventilation: Well-ventilated area
Eye Protection: Goggles
Protective Gloves: Yes - Latex
Protective Clothing: Avoid skin contact.

OHSA requires the use of eye protection & good ventilation when mixing, handling or discarding chemical disinfectants, & HICPAC Guidelines recommend the use of disposable plastic aprons & latex gloves when mixing, handling or disposing chemical agents.

Product Description:

Physical form: Translucent Liquid
pH: 12 Odour: Chemical

Date of Issue:
06/11/13



Enquiries to:

1 Bohnevista, 25 Payne Street, Pinetown P.O. Box 1338, Kloof, 3640, South Africa Tel: +27 31 702 7656
Fax: +27 31 702 3659 Website: www.steritech.co.za Email: info@steritech.co.za

