



SteriTech Disinfecting Instrument Precleaner

Concentrated Detergent-Disinfectant for cleaning instruments prior to sterilisation



SANS 636
Permit number
7333/10881



Available sizes:
1L; 5L; 20L
Dilution Factor:
1:20

Reg. No. ACT296NR529/242056/040/455

Quick:

Broad Spectrum disinfection in 5 minutes.

Duality:

Can be used for highly effective Manual Cleaning, and Ultrasonic Bath Cleaning.

Effective:

Contains a balanced system of non-ionic surfactants, sequestering agents, alkalis and grease cutters which removes both organic, and microbial loading. These prevent staining, and extend the lifespan of both the instruments, and the Cold Sterilant.

Easy:

Negates the need to scrub instruments in order to clean them, the soaking surfactant does all the work.

Protective:

Corrosion Inhibitors are added to protect instruments even when out of the liquid.

Antimicrobial Efficacy:

SteriTech Concentrate has been tested effective against the following organisms by the SABS, hospitals, government approved & internationally accredited laboratories in South Africa, UK and the USA:

Bactericidal: *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Escherichia coli* & Coliforms, *Chlamydia trachomatis*, *Salmonella typhi & choleraesuis*, *Enterobacter cloacae*, *Actinomyces naeslundii*, *Ureaplasma urealyticum*, *Streptococcus faecalis & mutans*, *Lactobacillus leichmanni & paracasei*, *Enterococcus hirae*, *Mycobacterium bovis & tuberculosis*, *Bacillus subtilis & cereus* and MRSA. **Virucidal:** Cytomegalovirus, HIV, Polio, Herpes, Influenza, Hepatitis, Rotavirus, Foot & Mouth Disease, Avian Influenza, SARS, Coronavirus, Norovirus, Parvo and Calicivirus. **Fungicidal** (including fungal spores): *Trichophyton mentagrophytes*, *Candida albicans* and *Aspergillus flavus & niger*.

Microbicidal: *Dermatophagoides pteronyssinus*



Active Ingredient:

SteriTech Disinfecting Instrument Precleaner, when diluted, contains SteriTech Concentrate (a Quaternary Ammonium Compound/TriButyl Tin Oxide Blend) at a level of 0.25%.

First Aid:

In case of eye contact:
Flush eyes thoroughly with water. Redness might occur, but should disappear within 1 hour. Should it persist, consult your physician.

In case of ingestion:
Induce vomiting, wash out mouth, then drink a large glass of milk mixed with a raw egg and consult your physician.

Instructions for Use:

Shake well before use. Dilute **SteriTech Disinfecting Instrument Precleaner** at 1:20 (eg 50mL product with 1L water). Instruments should be rinsed under running water & then thoroughly immersed in the solution, and soaked for 10-20 minutes. Agitation of the solution may be necessary if the instruments are heavily soiled. Remove instruments from solution and rinse thoroughly prior to sterilisation. Dilute only enough for immediate single use - if unused, discard within 7 days. Repeat as often as required by your Infection Control Protocol.

Recommended Uses:

- Cleaning and disinfection of Instruments prior to Sterilisation.
- Cleaning and disinfection of instruments possibly sensitive to scrubbing.
- Cleaning and disinfection of any metal equipment or surface which would benefit from the added Corrosion Inhibitor.
- Ultrasonic bath solution.

PPE Requirements:

Respiratory Protection: None
Ventilation: None
Eye Protection: None
Protective Gloves: Rubber or Neoprene for sensitive skin (for concentrate only).
Other Protective Clothing: None

Product Description:

Physical form: Colourless liquid
pH: 12-13 (conc.); 11-12 (dilute)
Rinsability: Excellent
Odour: Mild chemical odour
Incompatibility: anionic detergents
natural soaps
Stability: Stable



Enquiries to:

1 Bohnevista, 25 Payne Street, Pinetown P.O. Box 1338, Kloof, 3640, South Africa Tel: +27 31 702 7656
Fax: +27 31 702 3659 Website: www.steritech.co.za Email: info@steritech.co.za



SteriTech
RESIDUAL INFECTION PROTECTION