



SteriTech Disinfecting Suction Cleaner

Detergent-disinfectant for unblocking Suction Lines



SANS 636
Permit number
7333/10881



Available sizes:
1L; 5L; 20L
Dilution Factor:
1:100

Reg. No. ACT29GNR529/242056/040/496

Concentrated:

Economical to use in large quantities.

Residual Action:

Continues to disinfect even after the lines have been flushed, giving long-lasting protection, and ultimate peace of mind.

Effective:

Removes coagulated blood, saliva, organic tissue and other wastes to ensure that Suction equipment operates at peak efficiency at all times.

Faculty Tested:

Tested by both laboratories, and working facilities to ensure best formulation and results.



Antimicrobial Efficacy:

SteriTech Concentrate has been tested effective against the following organisms by the SABS, hospitals, government approved & internationally accredited laboratories in South Africa, UK and the USA:

Bactericidal: *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Escherichia coli* & Coliforms, *Chlamydia trachomatis*, *Salmonella typhi* & *choleraesius*, *Enterobacter cloacae*, *Actinomyces naeslundii*, *Ureaplasma urealyticum*, *Streptococcus faecalis* & *mutans*, *Lactobacillus leichmanni* & *paracasei*, *Enterococcus hirae*, *Mycobacterium bovis* & *tuberculosis*, *Bacillus subtilis* & *cereus* and MRSA. **Virucidal:** Cytomegalovirus, HIV, Polio, Herpes, Influenza, Hepatitis, Rotavirus, Foot & Mouth Disease, Avian Influenza, SARS, Coronavirus, Norovirus, Parvo and Calicivirus. **Fungicidal** (including fungal spores): *Trichophyton mentagrophytes*, *Candida albicans* and *Aspergillus flavus* & *niger*.

Microbicidal: *Dermatophagoides pteronyssinus*

Active Ingredient:

SteriTech Disinfecting Suction Cleaner, when diluted, contains SteriTech Concentrate (a Quaternary Ammonium Compound/TriButyl Tin Oxide Blend) at a level of 0.128%.

First Aid:

In case of eye contact:
Flush eyes thoroughly with water. Redness might occur, but should disappear within 2 hours. Should it persist, consult your physician.

In case of ingestion:
Do not induce vomiting, simply drink a large glass of milk mixed with a raw egg and consult your physician.

Instructions for Use:

Shake well before use. Dilute SteriTech Disinfecting Suction Cleaner at 1:100 (eg 5mL product with 500mL water). Pour or aspirate into the system. Allow a minimum of 5 minutes before flushing the system. For badly maintained systems (or the initial 4 uses), leave the solution in the system overnight, flushing the system the following morning. Dilute only enough for immediate use - if unused discard within 7 days of dilution. Repeat as often as required by your Infection Control Protocol (suggested Protocols available on request).

Recommended Uses:

- Disinfecting Suction equipment.
- Unblocking suction lines which contain congealed organic matter.
- Odour control in Suction Lines.
- Cleaning and Disinfection of hard to reach, non-sensitive areas in air conditioning systems.
- Can be used to unblock any tubing/piping where a build up of lime scale/calcium has occurred.

PPE Requirements:

Respiratory Protection: None
Ventilation: None
Eye Protection: None
Protective Gloves: Rubber or Neoprene.
Other Protective Clothing: None

Product Description:

Physical form: Clear Orange-Red liquid
pH: 13-14 (conc.); 11-12 (dilute)
Rinsability: excellent
Odour: Mild chemical odour
Incompatibility: anionic detergents
strong oxidisers
Stability: Stable

SABS
ISO 9001
SABS
ISO 13485

Enquiries to:

1 Bohnevista, 25 Payne Street, Pinetown P.O. Box 1338, Kloof, 3640, South Africa Tel: +27 31 702 7656
Fax: +27 31 702 3659 Website: www.steritech.co.za Email: info@steritech.co.za

SteriTech
RESIDUAL INFECTION PROTECTION