



SteriTech General Purpose Disinfectant

For disinfection of all pre-cleaned hygiene-sensitive areas



Environmentally friendly

Available sizes:
1L Trigger; 5L; 20L
Dilution Factor:
Ready-to-Use

Reg. No. ACT56NR529/242613/140/0672

Quick:

Bactericidal & virucidal in 1 minute.

Fungicidal in 2 minutes.


Effective:

Wide, practical range of tested organisms.

330 times more effective than traditional hypochlorite solutions.

Safe:

No typical 48hr half-life. Non-toxic.




Antimicrobial Efficacy:

SteriTech General Purpose Disinfectant contains Salvox at 1000ppm, which has been tested effective against the following standards & organisms:

Bactericidal: EN13727; EN1040; EN1276; SANS636; *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Escherichia coli*, *Streptococcus pyogenes*, *Enterococcus hirae*, *Bacteroides fragilis*, *Legionella spp* & *Listeria monocytogenes*. **Virucidal:** SANS6266; Polio virus & Adenovirus. **Fungicidal:** EN1275; EN1650; SANS636; *Candida albicans* and *Aspergillus brasiliensis (niger)*.

Active Ingredient:

SteriTech General Purpose Disinfectant contains Salvox (a stable, active hypochlorous acid solution) at a level of 1000ppm.



First Aid:

In case of skin contact: Wash immediately with plenty of soap and water.

In case of eye contact: Flush eyes with water for 15 minutes. Redness may occur, but should disappear within 30 minutes. Should it persist, consult your physician.

In case of ingestion: Do not induce vomiting, rinse out mouth, then drink a large glass of milk and consult your physician.

Instructions for Use:

Apply **SteriTech General Purpose Disinfectant** to a thoroughly pre-cleaned and rinsed surface, by spraying uniformly from a distance of approximately 30cm away. Ensure at least a 1 minute contact time for Bacteria & Viruses or a 2 minute contact time for Bacteria, Viruses & Fungi before the surface is re-used. Re-apply as often as required by your Infection Control Protocol - recommended protocols available on request.

Recommended Uses:

- Surface disinfection in hygiene-sensitive areas.

PPE Requirements:

Latex gloves & chemical splash goggles
But OSHA requires the use of eye protection & good ventilation when mixing, handling or discarding chemical disinfectants, and HICPAC Guidelines recommend the use of disposable plastic aprons & latex gloves when mixing, handling or disposing chemical agents.

Product Description:

Physical form: Translucent liquid
pH: 2.5-4.5

Odour: Faint chlorine

Incompatibility: Do not mix with other chemicals! Do not use on aluminium!

Date of Issue:
22/01/14



SteriTech
RESIDUAL INFECTION PROTECTION



Enquiries to:

1 Bohnvista, 25 Payne Street, Pinetown P.O. Box 1338, Kloof, 3640, South Africa Tel: +27 31 702 7656
Fax: +27 31 702 3659 Website: www.steritech.co.za Email: info@steritech.co.za

Please note: This is NOT a Material Safety Data Sheet