



SteriTech Prosthesis Gel

Disinfecting and Lubricating Gel for use on Prosthesis and Prosthetics



Available sizes:
250mL Push Disc; 5L
Dilution Factor:
Ready-to-use

Quick:

Disinfects on contact for instant infection control.

Limits friction:

Provides lubrication and so limits any friction between the Prosthesis and Prosthetic.

Soothing:

Decreases the discomfort created by stress, friction and infection on the amputation site.

Added Emollients:

Replaces lost natural skin oils to leave skin feeling soft and healthy.

Broad Spectrum:


Killing Bacteria, Viruses, Fungi & even fungal spores on contact, thus reducing the risk of infection to the amputation site, as well as eliminating any infection-based odours.

Antimicrobial Efficacy:

SteriTech Concentrate has been tested effective against the following organisms by the SABS, hospitals, government approved & internationally accredited laboratories in South Africa, UK and the USA:

Bactericidal: *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Escherichia coli* & Coliforms, *Chlamydia trachomatis*, *Salmonella typhi* & *choleraesius*, *Enterobacter cloacae*, *Actinomyces naeslundii*, *Ureaplasma urealyticum*, *Streptococcus faecalis* & *mutans*, *Lactobacillus leichmanni* & *paracasei*, *Enterococcus hirae*, *Mycobacterium bovis* & *tuberculosis*, *Bacillus subtilis* & *cereus* and MRSA. **Virucidal:** Cytomegalovirus, HIV, Polio, Herpes, Influenza, Hepatitis, Rotavirus, Foot & Mouth Disease, Avian Influenza, SARS, Coronavirus, Norovirus, Parvo and Calicivirus. **Fungicidal** (including fungal spores): *Trichophyton mentagrophytes*, *Candida albicans* and *Aspergillus flavus* & *niger*.

Microbicidal: *Dermatophagoides pteronyssinus*



Active Ingredient:

SteriTech Prosthesis Gel contains SteriTech Concentrate (a Quaternary Ammonium Compound/TriButyl Tin Oxide Blend) at a level of 1.25% & Isopropyl alcohol at a level of 60%.

First Aid:

In case of eye contact:
Flush eyes thoroughly with water. Redness might occur, but should disappear within 1 hour. Should it persist, consult your physician.

In case of ingestion:
Do not induce vomiting, wash out mouth, then drink a large glass of water and consult your physician.

Instructions for Use:

Apply a small amount of **SteriTech Prosthesis Gel** into the palm of your hand and massage onto the affected area. Do not rinse off. Also apply a small amount to the contact area of the prosthesis. Repeat application when necessary. Where fungal infection has already occurred repeat applications will assist in preventing cross contamination.
For external use only.

Recommended Uses:

- Lubricate the amputation site.
- Alleviate the discomfort created by friction.
- Prevent the occurrence of fungal infections.
- Deodorise the amputation site.
- Soothe friction-related redness & inflammation.

PPE Requirements:

Respiratory Protection: None
Ventilation: None
Eye Protection: None
Protective Gloves: None
Other Protective Clothing: None

Product Description:

Physical form: Clear Yellow Gel
pH: 6.5 - 7.5
Rinsability: excellent
Odour: Alcohol
Incompatibility: anionic detergents
natural soaps
Stability: Stable

SABS
ISO 9001
SABS
ISO 13485

Enquiries to:
1 Bohnevista, 25 Payne Street, Pinetown P.O. Box 1338, Kloof, 3640, South Africa Tel: +27 31 702 7656
Fax: +27 31 702 3659 Website: www.steritech.co.za Email: info@steritech.co.za



SteriTech
RESIDUAL INFECTION PROTECTION