

Vetguard Wound Protection Gel

A Liquid Plaster with Disinfecting Properties



Available sizes:
5mL Sachet; 250mL Push Disc
Dilution Factor:
Ready-to-use

Reg. No. G3766 (Act 36 of 1947)

Quick:

Dries to form a 'second skin' in 5-10 minutes depending on the type of wound.

Easy to Apply & Remove:

No distress to the wound bed & no training required to use it correctly.

Economical:

Saves on cost of plasters, bandages and dressings.

Residual Action:

Provides a Sterile dressing, and continues to disinfect even after skin has dried, meaning wounds stay contaminant free, giving long-lasting protection, and ultimate peace of mind.

Effective:

Keeps wounds moist, but contaminant free, resulting in accelerated healing, and decreased scar tissue formation.

Added Bitrex:

Discourages wound licking, and keeps flies away.



Vetguard Wound Gel

Reg No. G???? (Act 36 of 1947)

Wound Protection Gel for use
on Cats, Dogs & Horses

Store in a cool, dry place below 25 °C

Active Ingredients:
SteriTech Concentrate (Quarternary Ammonium
Compound / Tributyl Tin Oxide Blend) 1.25%
Gelling agents 2.50%
Inert ingredients 96.25%

Antimicrobial Efficacy:

SteriTech Concentrate has been tested effective against the following organisms by the SABS, hospitals, government approved & internationally accredited laboratories in South Africa, UK and USA: **Bactericidal:** *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Escherichia coli* & Coliforms, *Chlamydia trachomatis*, *Salmonella typhi* & *choleraesius*, *Enterobacter cloacae*, *Actinomyces naeshlundii*, *Ureaplasma urealyticum*, *Streptococcus faecalis* & *mutans*, *Lactobacillus leichmanni* & *paracasei*, *Enterococcus hirae*, *Mycobacterium bovis* & *tuberculosis*, *Bacillus subtilis* & *cereus*, MRSA and *Ornithobacterium rhinotracheale*.

Virucidal: Cytomegalovirus, HIV, Polio, Herpes, Influenza, Hepatitis, Rotavirus, Foot & Mouth Disease, Avian Influenza, SARS, Coronavirus, Norovirus, Parvo, Calicivirus, Infectious Bursal Disease virus and Newcastle Disease virus. **Fungicidal** (including fungal spores): *Trichophyton mentagrophytes*, *Candida albicans* and *Aspergillus flavus* & *niger*.

Microbiocidal: *Dermatophagoides pteronyssinus*.

Active Ingredient:

Vetguard Wound Protection Gel contains SteriTech Concentrate (a Quaternary Ammonium Compound/TriButyl Tin Oxide Blend) at a level of 1.25%.

First Aid:

In case of eye contact:
Flush eyes thoroughly with water. Redness might occur, but should disappear within 30 minutes. Should it persist, consult your physician.

In case of ingestion:
Do not induce vomiting, simply drink a large glass of water and consult your physician.

Instructions for Use:

Apply directly from dispenser to wound site. Spread into as thin a layer as possible using a gloved finger, gauze or spatula. Allow 5-10 minutes to dry into a skin-like layer. **DO NOT COVER THE TREATED AREA WITH A BANDAGE OR PLASTER.** To remove, wet area & wipe lightly with a clean cloth, cotton wool or gauze. Reapply daily or after washing or bathing. For external use only.

Recommended Uses:

- Use to cover wounds in hard to reach places or where scarring is unwanted.
- Excellent on the following wound types: Hot spots, bites, burns, debrided wounds, fragile skin, roasties, wounds in 'furry' places.
- Used as a liquid plaster/ second skin to keep wounds clean and bacteria free.

PPE Requirements:

Respiratory Protection: None
Ventilation: None
Eye Protection: None
Protective Gloves: None
Other Protective Clothing: None

Product Description:

Physical form: Clear Blue Paste
pH: 6 - 7
Rinsability: excellent
Odour: Mild chemical odour
Incompatibility: anionic detergents
natural soaps, strong oxidisers
Stability: Stable

SABS
ISO 9001
SABS
ISO 13485

Enquiries to:

1 Bohnvista, 25 Payne Street, Pinetown P.O. Box 1338, Kloof, 3640, South Africa Tel: +27 31 702 7656
Fax: +27 31 702 3659 Website: www.steritech.co.za Email: info@steritech.co.za



SteriTech
RESIDUAL INFECTION PROTECTION